ERAS APP for Patients



Final Presentation


Team Ramblin' on FHIR

Keith Adams, Lauren Aloia, Tyler Bobik, Chiawei Chien, Aaron Reed

Video Presentation Link: https://www.youtube.com/watch?v=j0mnaBZKMhE&feature=youtu.be

Background Research

- Bowel Surgery Handbook
- Sessions with mentors/physicians
- ERAS App study



An app for patient education and selfaudit within an enhanced recovery program for colorectal surgery: A pilot study on validity and usability.




₩ Tweet f 0 in 0

Nicolo Pecorelli, MD¹, Julio F Fiore Jr., MSc, PhD¹, Pepa Kaneva, MSc¹, Abarna Somasundram, BSc¹, Patrick Charlebois, MD², A S Liberman, MD², Barry L Stein, MD², Franco Carli, MD, MPhil², Liane S Feldman, MD². 'Steinberg-Bernstein Centre for Minimally Invasive Surgery and Innovation, ²McGill University Health Centre

Introduction: While patient engagement and clinical audit are key components of successful enhanced recovery programs (ERPs), they require substantial resource allocation. The objective of this study was to assess the validity and usability of a novel mobile device application for education and self-reporting of adherence for patients undergoing colorectal surgery within an established ERP.

Methods and procedures: Prospectively recruited patients undergoing colorectal surgery within an ERP received a tablet (iPad, Apple®) at their bedside during their hospital stay. Patients were instructed on how to use a novel app specifically designed to provide daily recovery milestones and record adherence to 15 different ERP processes and 6 patient reported outcomes (PROs) using a daily questionnaire for each day in hospital. Primary outcome was validity as measured by the agreement index (Cohen's kappa coefficient for categorical, and interclass correlation coefficient (ICC) for continuous variables) between patient reported data through the app and those recorded by a clinical auditor. Secondary outcomes included acceptability and usability of the app as measured by the System Usability Scale (SUS).

Results: Thirty-four patients participated in the study (mean age 61 years, 62% female, 62% laparoscopic surgery, 9% rectal surgery, 56% malignancy). Overall, patients completed 122 of 136 (90%) of the available questionnaires through the app. Median time to complete a questionnaire was 195 seconds (IQR 147–269). Substantial (kappa > 0.6) or almost-perfect agreement (kappa > 0.8) and strong correlation (ICC > 0.7) was found for 14 ERP processes and 5 PROs. Use of preoperative carbohydrate drink (kappa = 0.350) and ability to use the bathroom independently (kappa = 0.576) were the only elements with moderate agreement. Patient reported usability was high; median SUS score was 93 (IQR 83-99), with 79% of patients reporting an overall SUS score higher than 80 of 100. Scores were particularly high for the SUS learnability domain (median 100. IOR 81-100.) Only 5.



Innovation

Use of Ionic
 Framework
 allows cross
 platform
 capability



ios	Android	Windows
Setup: Name	Setup: Name	SETUP: NAME
Please enter your first and last name as shown on your medical record.	Please enter your first and last name as shown on your medical record.	Please enter your first and last name as shown on your medical record.
First Name	First Name	First Name
Last Name	Last Name	Last Name
Next >	NEXT →	Next >

Project Status

Working app that runs on iOS, Android and Windows phones.

• We were able to hit almost all the goals we planned on, this includes: Making a user friendly app that educates colorectal surgery patients, stores their information in a database and interacts with FHIR which I will explain more in depth in the demo.

Outstanding Elements

- Notifications that let the patient know if they miss a target.
- Save the profile as a questionnaire.


Deploy Final Product

• No actual deployment, we will just be handing the code off to our mentors, explaining to them how to set it up and run it. It is up to them if they want to put it up on an app store or not.


When you First Start the App Setup: Name

ios	Android	Windows
Setup: Name	Setup: Name	SETUP: NAME
Please enter your first and last name as shown on your medical record.	Please enter your first and last name as shown on your medical record.	Please enter your first and last name as shown on your medical record.
First Name	First Name	First Name
Last Name	Last Name	Last Name
Next >	NEXT →	Next >


Setup: Birthdate




Setup: Surgery Date



Setup: Smartphone




Setup: Support



Setup: Diabetic


GT FHIR Server Check for Diabetes on Marla Dixon Because it returns True, Yes is automatically checked.




If the user has no diabetes or the user's name is not located in the GT FHIR server the user must specify on their own whether they have diabetes or not.

Setup:	Diabetic
Are you Diabetic?	
Yes	
No	
< Back	Submit
31	

Main Menu and Education





Education: Pre-Surgery Screens





Pre-Surgery





Education: Day of Surgery Screens

Day of Surgery

Maintaining Body Temperature


A purple vest called "Bair Paws" will be


applied in the pre-op holding area. This

vest is meant to keep you warm prior to

surgery, and help maintain your


temperature during your procedure.
















Education: Day After Surgery Screens










Education: In-hospital Screens










Education: Recovery Timeline Screens





Recovery Timeline

(3 Month After Surgery) Typical

Recovery Timeline

You will likely have returned to work and

You will continue to get stronger and your appetite will improve

evaluate pouching due to changes in body shape and activity levels

Ostomy patients may need to re-

resumed regular activities







Targets: Pre-Surgery Targets

• Pre-Surgery targets are conditional on the surgery date and the current date


If <u>surgery date is more than 2 days before the current date</u>: Only these targets displayed **Pre-Surgery Targets**



Targets: Day Before Surgery

• Targets are conditional on the surgery date and the current date

If <u>surgery date is 1 day before the current date</u>: Only these targets displayed **Day Before Surgery Targets**



Time automatically added from when the user inputted it at the start via the database it was saved to.

Targets: Day of Surgery (Yes Diabetes)

Targets are conditional on two things:


- 1. The surgery date and the current date
- 2. If the patient has Diabetes or not: If Yes, then when the day of surgery matches the patient's **Day of Surgery** they entered on setup, drink a low-sugar sports drink card appears. If No diabetes is entered, the card will not appear.

If <u>surgery date = current date AND Diabetes = True</u>:

Only these targets displayed

Shows becuase user has Diabetes

Day of Surgery Targets



Targets: Day of Surgery (No Diabetes)

Targets are conditional on two things:


- 1. The surgery date and the current date
- 2. If the patient has Diabetes or not: If Yes, then when the day of surgery matches the patient's **Day of Surgery** they entered on setup, drink a low-sugar sports drink card appears. If No diabetes is entered, the card will not appear.

 If surgery date = current date AND Diabetes = False:

Only these targets displayed

ars because Day of Surgery Targets

Sports drink card Disappears because Diabetes is false



Targets: In-Hospital


• In-Hospital targets are conditional on the surgery date and the current date


If <u>surgery date was 1 day before the current date and less than a week after surgery</u>:


Only these targets displayed

In-Hospital Targets










Targets: Home Recovery

• Home Recovery targets are conditional on the surgery date and the current date


If <u>surgery date was 7 days before the current date</u>: Only these targets displayed **Home Recovery Targets**



Targets Integration With FHIR: Questionnair

For Questionnaire and QuestionnaireResponse on FHIR:


- 1. It saves when you leave targets.
- 2. Finds all the questions on the targets page that are 'radio'
- 3. Sends to server and checks if questionnaire already exists. If yes, it uses that questionnaire id.
- 4. Gets user's responses and creates QuestionnaireResponse resource. Sends to server and gets back id.




Targets Integration With FHIR: QuestionnaireResponse Resource

For Questionnaire and QuestionnaireResponse on FHIR:

- 1. it saves when you leave targets.
- 2. Finds all the questions on the targets page that are 'radio'
- 3. Sends to server and checks if questionnaire already exists. If yes, it uses that questionnaire id.
- 4. Gets user's responses and creates QuestionnaireResponse resource. Sends to server and gets back id.







QuestionnaireResponse Resource

```
"valueInteger":1
"valueInteger":1
"valueInteger":1
"valueInteger":2
```

Videos









Scroll


Profile







Contact












Gnatt Chart



Next Step (Future Improvement)

✓

- Deployment with client 12/7/2016
- Communication with Existing Client EHR Server
- Full SQLite Implementation
- Push notification from physician (ex: change of surgery date)
- Consult physician through in app messaging
- Admin login on physician end to track patient targets within app